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
MAINTENANCE
&
USERS MANUAL

THERMOBLOCK
DIGITALE

DIGITAL
THERMOBLOCK

FALC INSTRUMENTS s.r.l.

Treviglio (BG) Italy

INDICE

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INDEX

1.0 SECURITY.....	page 1
1.1 GENERAL.....	page 1
1.2 OPERATOR QUALIFICATIONS.....	page 1
1.3 INSTRUMENTS SECURITY.....	page 1
1.4 RESIDUAL RISKS.....	page 1
2.0 PRELIMINARY OPERATIONS.....	page 1
2.1 CHECK FOR DAMAGE DURING TRANSPORT.....	page 1
2.2 IN THE EVENT OF DAMAGE.....	page 1
2.3 CLEANING THE INSTRUMENT.....	page 1
3.0 INSTRUMENTS CHARACTERISTICS.....	page 2
4.0 COMMANDS.....	page 2
5.0 INSTALLATION AND USE PROCEDURES.....	page 3
5.1 PHYSICAL CHARACTERISTICS.....	page 3
5.2 FUNCTIONING.....	page 3
6.0 INTENDED AND UNINTENDED USES.....	page 4
7.0 MAINTENANCE.....	page 5
7.1 INSTRUMENTS MAINTENANCE.....	page 5
7.1.1 ORDINARY MAINTENANCE.....	page 5
7.1.2 PROGRAMMED MAINTENANCE.....	page 5
7.2 CLEANING.....	page 5
7.2.1 CLEANING PRODUCTS.....	page 5
7.3 TRANSPORT.....	page 6
7.4 STORING.....	page 6
7.4.1 STORING CHARACTERISTICS.....	page 6
8.0 PROBLEM RESOLUTIONS.....	page 6
9.0 ELECTRICAL WIRING DIAGRAM.....	page 7
10.0 DEMOLITION.....	page 10
10.1 OPERATOR QUALIFICATIONS.....	page 10
10.2 DEACTIVATING AND DEMOLITION.....	page 10
10.3 DEACTIVATING AND DEMOLITION PROCEDURES.....	page 10
10.4 RESIDUAL RISKS AFTER DEMOLITION.....	page 10

This manual should be conserved for the whole technical life of the instrument so that it can be easily consulted in case of need. If the instrument used should be sold, it should be sold complete of this manual.

Conserve the manual in a dry and easily reachable place for rapid consultation. Do not stain the manual to avoid creating stains that could prevent the lecture of this test.

1.0 SECURITY

1.1 GENERAL

Criteria and notions have been used in the construction of this instrument that guarantee that the instrument complies with the basic safety requirements under Low Voltage Directive 2006/95/EC and subsequent amendments. The accurate risk analysis carried out by the manufacturer means that the majority of risks connected with the use of the instrument have been eliminated. Full documentation of the adopted safety measures can be found in the technical file held by the manufacturer.

The manufacturer recommends that the instructions, procedures and recommendations contained in this manual, together with the security regulations in the working area, even for the use of the intended uses of protection.

FALC INSTRUMENTS cannot be held liable for any eventual damage to property, personal injury or injury to pets arising from the use of the instrument by unqualified operators and/or improper use of the instrument.

1.2 OPERATOR QUALIFICATIONS

This instrument must only be used by trained, qualified and authorised personnel fully familiar with the information contained in this manual.

1.3 INSTRUMENTS SECURITY

The instrument is fitted with a visual signal for the exceeding of the set point value. When the temperature is superior in respect to the set point of 2°C the alarm spy light  lights up and the heating is interrupted immediately.

This function is also used for reaching the set temperature.

1.4 RESIDUAL RISKS

Contact with the hot surfaces of the instrument such as the block holding tank, aluminium blocks or heated substances should be avoided by the users' limbs. We advise to use gloves or suitable tools for moving the hot containers.

ATTENTION

The non use or/and tampering with the protections could cause risks to the instruments users or other persons exposed; therefore *FALC INSTRUMENTS cannot be held liable for eventual damage to property, personal injury or injury to pets due to this behaviour*

2.0 PRELIMINARY OPERATIONS

2.1 CHECK FOR DAMAGE DURING TRANSPORT

Carefully remove the packaging and visually check the instrument for damage. Any dents in the visual parts indicate that the instrument has been knocked during transport, which could compromise its normal operation.

Check specially for the good conditions of the following parts:

- Instruments case.
- Controls and pilot lights.
- Check that all accessible screws and nuts are tight.

2.2 IN THE EVENT OF DAMAGE

Damage due to transport must be attributed to carrier and notified immediately to the manufacturer or his representative. Send a damage claim letter to the carrier.

2.3 CLEANING THE INSTRUMENT

Remove any dust or external dirt accumulated during transport.

3.0 INSTRUMENTS CHARACTERISTICS

Model		TD 150 P1	TD 150 P2	TD 150 P3
Places per block	n°	1	2	3
Maximum temperature	°C	+150	+150	+150
Precision per block	37°C	± 0,5	± 0,5	± 0,5
Heating power	W	90	180	270
External dimensions LxWxH	mm	220 x 185 x 95	220 x 225 x 95	220 x 325 x 95
Weight	Kg	2,2	2,8	3,6
Item		621.0422.01	621.0422.02	621.0422.03

ATTENTION → The maximum indicated temperatures refer to the internal part of block holder.


ATTENTION → The instrument's tare is available at 37°C

Power V/Hz 230/50-60	Security class 0	Protection class IP 40
-------------------------	---------------------	---------------------------

External casing in oven painted epoxy anti-acid powder paint, the block holder tank is in aluminium in order to assure a uniform heating distribution.

Control and setting of the values comes by means of a PID action microprocessor digital thermoregulator and PT100 probe for a better precision of the temperature. The heating of the thermoblock is set up by resistances underneath of the block holder tanks, all instruments are doted of a green light ON/OFF switch.

4.0 COMMANDS

For further information please see chapter
5.2 *FUNCTIONING* on page 3.

IG	Serves to switch on and off the instrument
Display	Visualizes the temperature value or the set temperature.
Enter	Confirms the modifications carried out.
Up	Increases the set temperature.
Down	Decreases the set temperature.
	Lights up when the resistances are in function or during the set phase.
	Lights up when the temperature is inferior to the set temperature.
	Lights up when the temperature is the same as the set temperature.
	Lights up when the temperature is superior to the set temperature.
	Lights up during the set phase or if the set point is over 2°C.

5.0 INSTALLATION AND OPERATING PROCEDURES

5.1 PHYSICAL CHARACTERISTICS

Position the instrument in a suitable place for use (norm D.L. 626/94-c risk valuation).

Other than the overall dimensions, the instrument also requires the following conditions:

- In the immediate vicinity of the instrument you must provide power sources. The instrument has only one electrical power point placed on the back of the same. Check that the electric power is suitable to the instruments power and that the same has a good ground earth (for tension information please see the label placed on the back of the instrument).
- Provide a sufficient flat space for the normal use and maintenance of the instrument, and any peripheral equipment

The instrument must be placed in a covered area, not exposed to the weather. Good lighting conditions for safe normal operation and maintenance of the instrument. The room where the instrument is used must be lighted sufficiently to allow for operation without risks due to areas in shadow. The accepted temperature for normal functioning is of a minimum of +5°C to a maximum of +40°C with maximum relative humidity of 80% for temperatures up to 31°C with linear reduction to 50% at 40°C.

ATTENTION

FALC INSTRUMENTS declines any responsibilities for inconveniences, malfunctioning and breakdowns which could be caused due to the non respect of the above mentioned indications.

5.2 FUNCTIONING

1) ELECTRICAL WIRING

The instrument has only one electrical power point placed on the back of the same.

2) TRIAL RUN

Before proceeding to the ordinary use of the instrument please control the perfect functioning of the same by carrying out one complete trial circle.

3) SWITCHING ON THE INSTRUMENT


- Place the power cable switch into the power point placed on the back of the instrument.
- Place the plug into the power socket.
- Place the general IG switch to position "1". The instrument will start heating.
-

4) SET PHASE

- By pressing Up or Down buttons the SET POINT value will appear.
- Use the Up or Down buttons to set the required temperature values.
- To memorise the new values press the Enter button.

If for four seconds no buttons are pressed the display will return to indicate the current temperature.

In the middle, underneath the display, there is a three LED movement indicator: green in the middle and red the two lateral ones. There are five possible signals:

	→ The temperature is inferior to 1,5°C from the SET POINT.
	→ The temperature is within -1°C from the SET POINT.
	→ The temperature is ± 0,5°C from the SET POINT.
	→ The temperature is within +1°C from the set point SET POINT.
	→ The temperature is higher of 1,5°C from the set point SET POINT.

5) SWITCHING OFF PROCEDURES

- Turn to position “0” the general IG switch.
- Remove from the instrument the containers and the blocks.
ATTENTION → movement of hot objects.
- Unplug instrument from the mains.

ATTENTION

FALC INSTRUMENTS cannot be held liable for eventual damage to property, personal injury or injury to pets arising from the non observance of the prescribed precautions or from the non use of the prescribed DPI's (See chapter *SECURITY* on page 1).

6.0 INTENDED AND UNINTENDED USES

😊 INTENDED USES

- 👍 The instruments are projected and built for dry heating and to maintain the temperature of solutions.
- 👍 The instruments are projected and built for treating substance normally used in laboratories.
- 👍 The instruments are projected and built to be used in laboratories for professional, industrial and teaching purposes.
- 👍 The instruments are usable in the biochemical, microbiological, haematology and diagnostics field.

😞 UNINTENDED USES

- 👎 Unintended uses are all those uses not clearly indicated in *INTENDED USES*.
- 👎 Use of the instruments in ambient conditions different to the indicated ones.
- 👎 Use of the instruments in explosive gas atmospheres.
- 👎 Electro medical uses.
- 👎 Washing the instrument with water jets.
- 👎 Use of acids or solvents in direct contact with the instrument.

ATTENTION

FALC INSTRUMENTS cannot be held liable for injury to property, personal injury, injury to pets arising from a use not in accordance with the above mentioned uses

7.0 MAINTENANCE

7.1 INSTRUMENTS MAINTENANCE

Carry out the maintenance operations with the instrument switched off and unplugged. Make sure that the instrument is cold before starting any type of operation.

Maintenance can be carried out in safety by expert personnel trained in and authorised for the use, preparation and maintenance of the instrument after having carefully read and fully understood the recommendations and instructions provided in this section.

7.1.1 ORDINARY MAINTENANCE

By ordinary maintenance we mean all operations that can be carried out by the user. These consist in cleaning operations, periodic and preventive inspections which allow the use of the instrument in complete safety.

7.1.2 PROGRAMMED MAINTENANCE

By programmed maintenance we mean all operation which can be carried out only by authorised personnel. These consist in periodic and preventive inspections on the instruments.

Nevertheless if the authorised and trained personnel should not be able to resolve a determined problem, contact the technical assistance.

7.2 CLEANING

ATTENTION

Carry out cleaning operations with the instrument switched off and unplugged from the mains.

After the end of each working day, thoroughly clean the instrument, especially the following parts:

- Plate and superior surface of the instrument.
- Frontal parts of commands.
- Surfaces surrounding the instrument.
- The floor around the working area.

7.2.1 CLEANING PRODUCTS

For cleaning operations it is sufficient to have a normal detergent and industrial water. Only use isopropyl alcohol in case of stubborn dirt.

ATTENTION

Avoid the use of solvents that can damage the synthetic material or the paint work. In particular avoid the use of benzene, perchlorate diluent and trichloroethylene.

7.3 TRANSPORT

If the instrument is transported by opened motor vehicles, the packaged instrument and its fittings must be covered completely by a plastic tarpaulin or equivalent to protect against the weather.

7.4 STORING

The indications contained in this section must be observed during temporary storage:

- ◆ Installing the instrument non immediately after being supplied..
- ◆ Deactivating the instrument and its storage whilst waiting for relocation.

7.4.1 STORING CHARACTERISTICS

- Up to 2000m. above sea level
- Accepted temperatures: from +5°C to +40°C.
- Maximum relative humidity 80% at temperatures up to 31°C with linear reduction to 50% at +40°C.
- Adequate weather-proofing.
- Sufficient room to manoeuvre for safe and easy lifting and transport.
- Horizontal surface capable of bearing more than the weight of the instrument.
- NEVER PLACE OR REST ANY OBJECT ON THE INSTRUMENT

8.0 PROBLEM RESOLUTIONS

Verified fault	Cause
The instrument does not work	4) The instrument is not connected properly at the mains. 5) The fuses placed on the back of the plug are burnt out. 6) There are some problems with the electrical line.
The instrument does not heat	5) The instrument is switched off. 6) The value of the set temperature is inferior to the ambient temperature. 7) * The instruments digital thermo regulator does not work properly (contact technical assistance). 8) * The resistances are burnt (contact technical assistance).


**** The above described faults are attributed to the use of the instrument without using the manual procedures, to the technical life of the various components and long time use.***

In the event that the authorised and trained personnel are not able to resolve a determined problem please contact the technical assistance.

9.0 ELECTRICAL WIRING DIAGRAMM


MODEL : TD 150 P1

NR	COMPONENTS	NR	COMPONENTS
A	FUSES	D	TM FALC THERMOREGULATOR
B	COUPLING SOCKET	E	PT100 PROBE
C	TWO POLE SWITCH	F	RESISTANCE


MODEL : TD 150 P2

NR	COMPONENTS	NR	COMPONENTI
A	FUSES	D	TM FALC THERMOREGULATOR
B	COUPLING SOCKET	E	PT100 PROBE
C	TWO POLE SWITCH	F	RESISTANCE

MODEL : TD 150 P3

NR	COMPONENTS	NR	COMPONENTS
A	FUSES	D	TM FALC THERMOREGULATOR
B	COUPLING SOCKET	E	PT100 PROBE
C	TWO POLE SWITCH	F	RESISTANCE

10.0 DEMOLITION

10.1 OPERATOR QUALIFICATIONS

A qualified operator who has read and fully understood the contents of chapter *SECURITY* on page 1 and the contents of chapter *DEMOLITION* on page 10.

10.2 DEACTIVATING AND DEMOLITION

Once the instrument has reached the end of its technical operation life, the instrument has to be deactivated. Putting the same out of function and putting the instrument in a condition not to be used any more for the scopes which it had initially been projected and constructed, it should, however, be possible for the raw materials of the same to be re-used.

It is necessary for a deactivating and demolition to be carried out in security of the instrument.

ATTENTION

FALC INSTRUMENTS cannot be held liable for personal injury or injury to pets arising from the re-use of the single parts of the instrument or installation situations different to the original ones.

10.3 DEACTIVATING AND DEMOLITION PROCEDURES

The activating and demolition procedures of the instrument have to be carried out only by fully qualified and trained personnel.

Switch off the instrument.

Disconnect the power sources.

Dismantle the following parts:

- Electrical parts.
- Plastic parts.
- Resin Support Feet.

*In case of movement of the instrument refer to chapter *INSTALLATION AND OPERATING PROCEDURES* on page 4, *TRANSPORT* and *STORING* on page 7.*

10.4 RESIDUAL RISKS AFTER DEMOLITION

If the indications of section *DEACTIVATING AND DEMOLITION* page 10, are accurately carried out, there are no residual risks after demolition.

ATTENTION

The instrument has been constructed with non biodegradable material. Dispose of the instrument in an authorised deposit for the disposal.

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FALC INSTRUMENTS s.r.l.

SUBJECT: GUARANTY

For a period of 12 months from the date of delivery to the carrier. With the word "guaranty" we engage ourselves to substitute any and every faulty instrument, product components or part, due to faults or defect of our manufacturing process, not due to transport or misuse, as soon as we are informed and the item sent back to us for verification. The item should be sent back to us at your charge with a report explaining and describing the fault. When the problem has been verified and confirmed, the item shall be substituted or changed repaired and sent back to the customer free of charge. In a case the item has been tampered with, misused, or used not following the instructions, or in case wrong repair trials have been done on it, there shall be no guaranty whatsoever.

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The following unit complies with the safety requirements outlined by the EC DIRECTIVE concerning the guidelines for ELECTROMAGNETIC COMPATIBILITY (2004/108/EC) and LOW VOLTAGE REGULATIONS (2006/95/EC).

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The units are manufactured in compliance with the regulations

CEI EN 61000-6-2

CEI EN 61000-6-4

CEI EN 61010-1

CEI EN 60355-1

CEI EN 61326-1

Il responsabile legale / legal responsible

Scotti Francesco